

COLORIMETRIC / SPECTROPHOTOMETRIC PROCEDURES

Colorimetric/Spectrophotometric recording requirements in the **Sample Analysis Records** should include:

- a) wavelength (if adjustable)
- b) cell path length (if interchangeable)
- c) lot number of commercially prepared reagent in use
- d) calibration curve referenced (or slope, y-intercept, and correlation coefficient)
- e) date, time and analyst (sample analysis)
- f) some methods require logging date and time sample collected due to short holding time (ex. chlorine)
- g) instrument identification (if more than one is available)
- h) type of parameter/analyte
- i) some methods may require temperature
- j) program area (drinking water, wastewater, or SHW) - when lab performs parameter analysis for more than one Act
- k) sample identification
- l) sample results (with units of measurement: ex.: mg/L, etc). Some methods may require calculations to attain sample results. If so, the laboratory will need to record sample volume, absorbance, and any other factor used in the calculation.
- m) comment section

Volume, concentration of spike, spike reading, and percent recovery must also be recorded in the **Sample Analysis Records** (when applicable). In most methods (excluding chlorine) spikes and duplicates must be performed at a frequency of 10%.

Daily calibrations of the instrument must include, at minimum, a blank and two standards.

The **calibration curve** must be established per lot of reagent(s) for each instrument used (ex.: if free and total residual determinations are made, a curve must be prepared for each, due to different lots of DPD reagents in use), either annually or when a new lot of reagents are used. Some methods may require more frequent curves. This curve must be comprised of at least three (3) standards and a blank that cover the range of the method although five to seven standards are recommended (refer to each method, this requirement may vary). A laboratory that is using linear regression to define a curve must use at least five (5) standards. These values of true concentrations versus observed readings, or absorbance versus true concentration must be documented with the curve.

The **calibration curve** should include:

- a) date generated
- b) analyst's initials
- c) instrument identification
- d) lot number of reagent/standard
- e) cell path length
- f) wavelength
- g) true values of standards and blank versus observed values of standards and blank (or absorbance versus concentration)
- h) graph of true values (x-axis) and observed values (y- axis). Some methods/parameters may need to use absorbance versus concentration and may allow for the "physical" (plotted) curve to be substituted for the following logged information: y-intercept, slope, and correlation coefficient
- i) parameter method number and reference (SM, EPA, etc.)

This curve should be a linear regression line (or "best fit" line) and analysts should not "connect-the-dots".

For those samples requiring pretreatment (digestion, distillation), a blank and at least one standard should also be pretreated. Sample dilutions must be made prior to color development (an additional sample portion will be needed for this purpose).